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CERVICAL AND DORSAL HEATING PAD

PRODUCT DESCRIPTION

1. Heating surface
2. Lower adjustment strap
3. Upper adjustment strap
4. Detachable temperature control
5. Temperature selector
6. Temperature display

SAFETY INSTRUCTIONS

Carefully read the instruction manual.

Protection category II.

Do not pierce the heating pad!

Do not use the heating pad when it is folded!

It must not be used by very young children (0-3 years)!

Only use the heating pad indoors!

Maximum washing temperature
30 °C, very gently wash.

Do not chlorinate!

Do not dry the heating pad in the dryer!

Do not iron the heating pad!

Do not dry clean!

- **Please read this instruction booklet carefully before starting the appliance and store it for subsequent consultation. Non-observance and non-compliance with these instructions could result in an accident.**

- Carefully check the heating pad before each use for any signs of wear or tear.

- Do not use if there are any signs of damage or improper use in the appliance, the switch or the cables, or if the device does not work.

- Before connecting the device to its power source, make sure that the power voltage indicated on the classification plate is compatible with your power supply.

- Do not fold or wrinkle the heating pad during use.

- Do not pierce the pad, or place safety pins or other sharp objects in the pad.

- Do not use the device on children or on people who are in any way challenged, asleep or insensitive to heat (people who

are unable to react to overheating).

- Children under the age of 3 years cannot use this device because they are unable to react to overheating.
- The device must not be used by children over the age of 3 years, unless the control unit has been adjusted by parents or legal tutors, or unless the child has been clearly instructed on how to use the control unit correctly.
- This device can be used by children aged 8 and over and by people with reduced physical, sensory, or mental abilities, or those lacking experience and knowledge, provided they have been given supervision or instructions on how to use the device safely and understand the risks involved.
- Children shall not play with the appliance.
- Cleaning and maintenance of the device must not be carried out by unsupervised children.
- Never use the device without supervision.
- The device is only for private use and must not be used in hospitals.

- Do not fall asleep while the heating pad is switched on.
- Prolonged use of the heating pad at a high setting could cause burns on the skin.
- Do not cover the pad with another cushion. The control switch must not be covered or placed on or under pad while the unit is working.
- Do not sit on the heating pad.
- Never touch a heating pad that has fallen into water. Unplug the unit from the outlet immediately.
- Keep the power cable away from hot surfaces.
- Never carry, pull or turn the heating pad by the power cable and make sure the cable does not get tangled.
- Do not use the heating pad when it is wet and only use it in dry environments (not in bathrooms and similar environments).
- The controls and cables must not be exposed to any kind of humidity.
- If there is a fault, do not try to repair

the unit yourself. Repairs should only be carried out by an authorised specialised electrician* or other appropriately qualified personnel.

- If the power cable is damaged, it must be replaced by the manufacturer, or similarly qualified persons to prevent any danger.
- A damaged power cable must be replaced by the same type of power cable.
- You are recommended to check the appliance for damage from time to time. Never use the appliance if the cord or the appliance show signs of damage. If such is the case, repairs must be carried out by a qualified electrician (*). If the cord is damaged, it must be replaced by a qualified electrician (*) to avert any danger.
- The heating pad can only be operated with the control unit provided, and not with any other.
- Leave the heating pad to cool before folding and storing it.
- If you have any concerns over health

Competent qualified electrician: Any person who is qualified, authorised and competent to do this type of repairs with an aim to avoid danger of any kind.

problems, consult your doctor before using the heating pad.

- Do not use the heating pad on any part of the body if it is swollen, injured or has any inflammation.
- If you have persistent pain in muscles and joints, let your doctor know. Persistent pain could be a system of serious illness.
- If you feel pain or unease while using the unit, stop using it immediately.
- This appliance is not a medical device, which is why it does not for health care purposes.

HOW TO USE

- Remove all the packaging from the appliance and make sure that it is complete and that no damage has been caused during its transport.
- The appliance comes with an unattached cable. Before use, connect the cable to the pad by exerting slight pressure and the other end of the plug to a power outlet.
- The control switch has 3 temperatures. Select the required one by sliding the selector in one of the following positions

1_minimum temperature ($53^{\circ}\text{C} \pm 5^{\circ}\text{C}$)

2_moderate temperature ($59^{\circ}\text{C} \pm 5^{\circ}\text{C}$)

3_maximum temperature ($65^{\circ}\text{C} \pm 5^{\circ}\text{C}$)

- The thermal sensation will depend on different factors such as the ambient temperature, the clothes you are wearing, the pressure of pad during its use, etc,...
- For quick heating, select the maximum temperature level and then set the temperature to your requirements.
- The device will switch off automatically 90 minutes after it is switched on, as a safety measure.

CLEANING AND STORING

- Before cleaning the pad, unplug it from the mains and leave it to cool down for at least ten minutes.
- The heating pad is fitted with a detachable cable. Unplug the plug and remove the cable from the heating pad.
- You can wash the heating pad in the washing machine, maximum 30°C and with a gentle spin.
- Keep the unit in a clean and dry place, fully extended and with nothing covering or on top of it.
- Only use the appliance once it is completely dry.

TECHNICAL SPECIFICATIONS

Model: 5427054

Size: 63x42cm

Rated voltage: 220-240V ~ 50-60Hz

Power: 100W

This appliance conforms to Directive 2014/35/EU on Low Voltage and with Directive 2014/30/EU on Electromagnetic Compatibility.

ENVIRONMENTAL PROTECTION - / EU DIRECTIVE 2012/19

In order to preserve our environment and protect our health, any waste electrical and electronic equipment must be disposed of in accordance with specific rules with the implication of suppliers and users.

The symbol on the packaging indicates that this product must not be disposed of together with ordinary household waste. The user must take the product to a collection point for the recycling of electrical and electronic appliances.

Disposal must be in line with prevailing local regulations. Consult your regional or municipal authorities about the options for waste disposal.

GUARANTEE TERMS

The legal guarantee conforming to prevailing legislation is applicable to this product. For the duration of the legal guarantee, all material and manufacturing defects are covered by the guarantee.

To avail of the guarantee, you will need to present the purchase receipt at the store where the item was bought. You can then choose either a refund or an exchange for a new item.

For after-sales service contact us with any query you may have at: <http://adibera.com/es/contact/questions>

No claims will be accepted for broken accessories or defects

due to inappropriate use or care, excessive wear due to intensive, non-domestic use, incorrect use due to failure to comply with the instruction manual, or normal wear. If this appliance is used incorrectly, for professional or semi-professional purposes or in a way that is not in line with the operating instructions, the guarantee will be null and void.

No guarantee claims will be accepted if the appliance has been modified or if any attempts have been made to repair it outside the official technical assistance services. In such an event, all liability in the event of accidents will be rejected.

The parts subjected to wear, as well as cleaning and maintenance or the replacement of worn parts are not covered by the guarantee and must be paid for by the buyer.

This instruction manual and any future updates can be digitally downloaded at www.adibera.com

AFTER THE LEGAL PERIOD OF GUARANTEE

Repairs after the guarantee period will be on the owner's account.

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